

## Piezo Mini

# MINI DENTAL PIEZO SURGERY

Minimize and lightweight  
dual-screen design for smart dual control



# MINI DESIGN

## COMPACT & EFFICIENT

# DentalDeal

Balancing minimally invasive techniques & power output

## 5"

MAIN SCREEN



Full-touch interface for intuitive operation

## 1.6"

ROUND SCREEN



Real-time data display and adjust

### OPTIMIZED DIMENSIONS

perfectly designed for a streamlined workspace

### SLIM HANDPIECE

Operates at full power for 10 mins with a handpiece temperature of only 33°C, ensuring safety and comfort



Bone Mode



Endo Mode



Perio Mode



Custom Mode



Clean Mode

### PARAMETER

**VOLTAGE** 100~240V 50/60Hz 120VA

**FUSES** 2×F10A 250V

**FREQUENCY** 26~38Hz

**RIDAL AMPLITUDE** 20-200µm

**LATERAL AMPLITUDE** ≤30µm

**WATER FLOW** 60-100ml

[www.dentaldeal.ca](http://www.dentaldeal.ca)

PLUWA<sup>®</sup>  
PLUS POWER  
SG60P  
SN: SU2524JA11823