

## R1 Plus Gutta-percha Cutter



# DentalDeal


Designed by international dental experts

Conducting 220 precise manufacturing processes on the equipment imported from Germany and Japan




## 01

Gutta-percha point trimming,  
smooth and glabrous cut surface



## 02

18 trimming and diameter measuring  
holes with an accuracy of 0.01mm,  
suitable for all gutta-percha points  
on the market



## 03

Anti-rust and double-sided cutting  
blade, smooth cutting can be done  
by sliding forward or backward



## 04

30mm length measuring ruler



## 05

Laser marking guarantees durable  
and fadeless marks



## 06


Made of aluminum alloy,  
autoclavable, light weight, no rust



Guilin Woodpecker Medical Instrument Co., Ltd.  
Information Industrial Park, National High-tech  
Zone, Guilin, Guangxi, 541004 P. R. China.

Registrant/manufacturer/after-sales service provider:  
Guilin Woodpecker Medical Instrument Co., Ltd.  
Sales Dept.: 0773-5873196/2350599  
Fax: 0773-5822450  
After-sales Service Dept: 0773-5827898/13978361362  
E-mail: woodpecker4@glwoodpecker.com  
Website: <http://www.glwoodpecker.com>

Dealers



ZMN-XC-035